

Quick Start Guide

The LaserComb is a very effective product if used correctly. It is manufactured using precision components and should last for many years under normal usage.


We recommend using the LaserComb for a minimum of 5 minutes or a maximum of 15 minutes per treatment.

- 3 times a week for the first 4 weeks
- 2 times a week thereafter

It is recommended but not essential that the LaserComb be used on a clean scalp after washing. Your hair can be wet, towel dried or dry. The LaserComb should be used before applying any other topical treatments, gels, sprays or concealing products.



1. Insert the power attachment into the LaserComb.



2. Plug in the wall adaptor.



3. Turn on the LaserComb by first switching on the 'Power' button followed by the 'Laser' button.



4. Place the LaserComb flat on your scalp using the two rows of teeth.



Stop and holding 4-5 seconds every 1/2 inch.

5. Move the LaserComb using a 'Step' method by placing the LaserComb on a spot and leaving it there for 4 or 5 seconds then moving it 1/2 inch to the next spot. (Mentally counting "one thousand, two thousand, etc." should help with timing).

- 6 Move the LaserComb from the front of your hairline to the back, being sure to cover the whole scalp. Then from bottom to top, again covering the whole scalp. When finished, you will have covered your scalp twice. With long hair, use a brush or your hand to lift your hair in advanced of the LaserComb so the laser light reaches the scalp more effectively.
- 7 After treatment, turn off the power button, remove the wall plug and return the LaserComb and power supply to the case for safekeeping.

Please Note:

- The total time it should take will range between 5 and 15 minutes.
- It is important to move the LaserComb very slowly. If you move the LaserComb too fast, the treatment will not be as effective.
- It is recommended that you do not overuse the LaserComb. As with any therapy, as there is a point of maximum return where the best results are achieved.

Please refer to the user manual for additional product information.